TempoCemID

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English

Product description

TempoCemID is a translucent, dual-cure composite luting cement that is indicated for temporary luting of temporary crowns, bridges, veneers, inlays and onlays, as well as all kinds of temporaries. The material is available in the Smartmix-syringe.

Indications

- Temporary luting of temporary crowns, bridges, inlays and onlays.
- Temporary luting of temporary veneers.
- Temporary luting of implant-borne temporary prosthesis.
- Semi-permanent luting of temporary crowns, bridges, inlays, and onlays.
- Semi-permanent luting of implant-borne temporary prosthesis.

Contraindications

Do not use the material in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

Basic safety information

- For dental use only.
- Keep out of the reach of children.
- Avoid contact with the skin. In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with the eyes. In the event of accidental contact with the eyes, immediately rinse eyes thoroughly with plenty of water and consult a physician if necessary.
- In the event of application to dentine in close proximity to pulp, use suitable pulp protection.

Side effects

There are no known systemic side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

Interactions

Products containing eugenol, for example temporary cements or temporary root canal filling materials, may have an effect on the curing of the material. Products containing eugenol and other substances that inhibit polymerization must not be used together with the material.

Timing

Working time (including mixing time)*	1:00 min
Removal of excess material (without light-curing)	1:15–2:00 min after insertion of the restoration
Light-curing time for removal of excess (optional)	1–2 s per side
Intraoral curing	6:00–7:00 min
Light-curing	20 s per side

*The specified times apply at a room temperature of 23 °C/73 °F and a normal relative humidity of 50%. Higher temperatures shorten these times and lower temperatures increase them.

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material. Observe the manufacturer's instructions for the light-curing unit.
- The material sticks to composites. Apply a thin layer of separating agent (e.g. vaseline, K-Y Jelly) to existing composite-based fillings, underfillings or build-up fillings before using the material. In the case of further use of the temporary, apply separating agent in advance.
- When using a syringe for the first time, express a pea-sized amount of material and discard it.
- The removal of the cement from the restoration is made easier by an ultrasound.

Recommended use

1. As required, apply separating agent to the tooth core, core build-up, or temporary (see also Notes for use).
2. Apply the material to the temporary.
3. Set the restoration within 1:00 min after mixing has begun and attach to the dry, prepared tooth using slight pressure.
4. If necessary, instruct the patient to adopt habitual occlusion.

- Remove excess cement within 1:15 to 2:00 min of seating the restoration with the aid of a microbrush, brush, foam pellet, or scaler. Drying with oil-free and water-free air makes excess cement easier to identify. Remove interproximal excess by flossing carefully. It is alternatively possible to remove excess cement by briefly curing it using light (1–2 s per side) and then detaching it with a scaler.
- Allow the material to chemically cure completely for 6:00 to 7:00 min after mixing has begun.
- Translucent restorations: Additionally cure the material with a suitable light-curing unit for at least 20 s per side, as required.

Use of the Smartmix-syringe

- Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
- Attach a new mixing tip [B].

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix-syringe [C].

- Stop the mixing tip by turning it 90° clockwise.
- The material is mixed in the cannula when it is expressed and it can be applied directly.

Note: When using a syringe for the first time, express a pea-sized amount of material and discard it. After application, leave the used mixing tip on the Smartmix syringe as a cap.

Composition

Cement: Glass filling material in a matrix of multifunctional methacrylates. Catalyst, stabilizer, additives. Does not contain methyl methacrylate or peroxides. Filler degree: 40 wt % = 22 vol.% (0.02–2.5 µm)

Storage and shelf life

- ▶ Store at room temperature (15–25°C/59–77°F) in a dry place.
- ▶ Keep away from direct sunlight.
- ▶ Do not use after the expiration date.

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